

OUTPUT CHARACTERISTICS

Maximum pressure:.....90 mm Hg
Number of cameras:.....24
Treatment time:.....0 - 60 min
Retention time:.....2 - 4 - 6 sec
Rest time:.....4 - 6 - 8 sec
Classification:..... Class I type BF

Protection fuses:.....2 x T1,5A / 250 V
Degree of protection:..... IPX0
Preset programs:.....6
Customizable programs:..... Yes
Autonomous treatments:..... Yes

TECHNICAL CHARACTERISTICS

Supply voltage:.....120 - 230 V 50/60 Hz
Maximum consumption:.....180 VA (0.75 A at 230 V)
Interface:.....color screen 10.4 "Touch screen

EXCLUSIVE ITALIAN DISTRIBUTOR



SPORTARREDO GROUP S.C.
Via dell'Industria, 7 - 30020 GRUARO (VE)
Tel. +39 0421 767676 - Fax. +39 0421 767670
www.sportarredo.com - info@sportarredo.it

Sapphire

KiesPress²⁴

PRESSOMASSAGE
REDUCE, MOLD AND ELIMINATE.



With guarantee



Sapphire

KiesPress

Stimulates the circulatory system

The pressomassage is a treatment that uses the air pressure to perform a lymphatic drainage.

Generates Lymphatic Drainage

With Sapphire KiesPress24, oxygenation improves throughout the body, obtaining a great benefit for the entire body.

Counteract cellulite, edema and lymphedema

Correcting of alterations in the circulatory system and imperfections of cellulite, reducing it and preventing future training.



Reduce, Mold and Eliminate

With Kies Press24, you will stimulate the circulatory system, generating a lymphatic drainage and contributing to the elimination of cellulite, edema, lymphedema and venous disorders such as varicose veins.

The pressomassage is a treatment that uses the air pressure to perform a lymphatic drainage. It improves oxygenation throughout the body. This effect being a great benefit for the whole organism.

Perhaps the best place to check is the skin, whose appearances will be improved in terms of tone and firmness. Moreover, increasing the oxygen in the blood also increases the tissue nutrition and facilitates cell regeneration.

EQUIPMENT DELICATELY DESIGNED AND MADE JUST FOR YOU

Made in Spain

KiesPress24 device is capable of achieving a maximum level of your needs and perfect results.

PROGRAMMABLE

Digitally regulated
WITH 6 OPTIONS OF TREATMENT

TYPES OF MESSAGES

ALL PROGRAMS HAVE 2 MESSAGE OPTIONS:

WAVE EFFECT
SEQUENTIAL

It has a digitally regulated monitor, where all the great functions of this new hi-tech pressotherapy systems are found.

The device has 6 programs to choose from and are ideal for all types of sessions as well as two modalities of massage: Wave and Sequential.

With guarantee
Sapphire

